

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Polyamid 6 Pulver reinst	Kat.Nr./Cat.No. :	33143
Nylon 6, research grade	Lot/Contr.No. :	180575

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weißes Pulver white powder	entspricht corresponds
Korngröße (mm) particle size		0.05 – 0.16	entspricht corresponds
pH	10 % in H ₂ O	6.0 – 7.0	6.3
Mindesthaltbarkeit minimum shelf life			23.05.2026
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 14.05.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.