

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Ponceau S Lösung für Elektrophorese (0.2%)	Kat.Nr./Cat.No. : 33427
Ponseau S solution for electrophoresis (0.2%)	Lot/Contr.No. : P130076

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		rote Lösung red solution	entspricht corresponds
Konzentration concentration		0.2 % in 3 % TCA	entspricht corresponds
A 1 cm/ 524 nm	1:100 verd. in H ₂ O 1:100 diluted in H ₂ O	min. 0.9	1.034
Proteinfärbung protein staining	CAF	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			03/ 2016
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 07.03.2013

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.