

CERTIFICATE OF ANALYSIS

Proteinase K from Tritirachium album lyophil.	Cat.No. : 33752 Contr.No.: 130986
---	--

Parameter	Method	Specification	Result
Molecular weight		ca. 28 000	
Appearance		white lyophilisate	corresponds
Activity	DMC U/mg determined with TNBS	min. 8	13.3
Extraneous activities	DNase RNase	not detectable not detectable	corresponds corresponds
Minimum shelf life			09/ 2015
Storage (°C)			+2 to +8

Unit definition

One DMC-Unit is that amount of enzymatic activity which catalyzes the cleavage of 1 μ equivalent peptide bond from dimethylcasein per minute at 25°C, pH 7.0, expressed in terms of the appearance of new terminal amino groups.

We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.
--

SERVA Electrophoresis GmbH
Quality Control**Printing date:** 04.09.2013

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.