

CERTIFICATE OF ANALYSIS

Proteinase K from Tritirachium album lyophil.	Cat.No. : 33752 Contr.No.: 180916
---	--

Parameter	Method	Specification	Result
Molecular weight		ca. 28 000	
Appearance		white lyophilisate	corresponds
Activity (mAnsonU/mg)		min. 30	≥ 40
Extraneous activities	DNase	not detectable	corresponds
	RNase	not detectable	corresponds
Minimum shelf life			30.09.2020
Storage (°C)			+2 to +8

Unit definition

One unit is defined as the amount of enzyme that liberates Folin-positive amino acids and peptides, corresponding to 1 pmole tyrosine under assay conditions in 1 minute using hemoglobin as substrate.

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 09.09.2019

Christian Monsler

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.