

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS**SERVA Violet 17 Staining Kit****Kat.Nr./Cat.No. : 35074****Lot/Contr.No. : P240376**

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Anwendung</b> application	zum Färben von IEF-Gelen for staining of IEF gels	entspricht corresponds
<b>Inhalt</b> content	1 g SERVA Violet 17 500 ml 20% TCA Solution for Fixation 500 ml Stain Solubilizer Solution (20% phosphoric acid) 2 x 500 ml Destaining Solution (3% phosphoric acid)	
<b>Mindesthaltbarkeit</b> minimum shelf life		10.03.2026
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**

**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control****Druckdatum/ printing date: 02.12.2024**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**

This report has been computer generated and does not contain a signature.