

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

| | |
|---|--|
| SERVA Silver Stain Kit for SDS PAGE | Kat.Nr./Cat.No. : 35076 Lot/Contr.No. : P250018 |
|---|--|

| Parameter parameter | Spezifikation specification | Ergebnis result |
|---|---|---------------------------|
| Anwendung application | SDS PAGE | entspricht corresponds |
| Inhalt content | 500 ml Fixing solution 25 x 30mg Sodiumthiosulfate 5 H ₂ O 250 ml Silver nitrate solution 1.5 ml Formaldehyde 500 ml Sodium carbonate solution 30 g Glycine | |
| Mindesthaltbarkeit minimum shelf life | | 31.03.2026 |
| Lagerung (°C) storage | | +15 bis/ to +30 |

**Bei Lagerung unter +10 °C kann die Natriumcarbonat Lösung auskristallisieren.
Die Kristalle lösen sich beim Erwärmen des Produktes wieder.**

Storage below +10 °C may cause crystallization of Sodium carbonate solution.
Crystals can be solved by warming up the product.

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 23.01.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.