

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS**SERVA Silver Stain Kit**  
for Native PAGE**Kat.Nr./Cat.No. : 35077**  
**Lot/Contr.No. : P130308**

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Anwendung</b> application	native PAGE	
<b>Inhalt</b> content	500 ml Trichloroacetic acid 25 x 30mg Sodiumthiosulfate 5 H <sub>2</sub> O 250 ml Silver nitrate solution 1.5 ml Formaldehyde 500 ml Sodium carbonate solution 30 g Glycine	
<b>Mindesthaltbarkeit</b> minimum shelf life		05/ 2015
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control****Druckdatum/ printing date: 03.03.2014**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.