

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS**SERVA Silver Stain Kit**
for Native PAGE**Kat.Nr./Cat.No. : 35077**
Lot/Contr.No. : P150298

| Parameter parameter | Spezifikation specification | Ergebnis result |
|---|---|---------------------------|
| Anwendung application | native PAGE | |
| Inhalt content | 500 ml Trichloroacetic acid 25 x 30mg Sodiumthiosulfate· 5 H ₂ O 250 ml Silver nitrate solution 1.5 ml Formaldehyde 500 ml Sodium carbonate solution 30 g Glycine | |
| Mindesthaltbarkeit minimum shelf life | | 22.06.2016 |
| Lagerung (°C) storage | | +15 bis/ to +30 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control**Druckdatum/ printing date: 28.08.2015**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.