

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

SERVA Silver Stain Kit
for Native PAGE

Kat.Nr./Cat.No. : 35077
Lot/Contr.No. : P150374

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	native PAGE	
Inhalt content	500 ml Trichloroacetic acid 25 x 30mg Sodiumthiosulfate· 5 H ₂ O 250 ml Silver nitrate solution 1.5 ml Formaldehyde 500 ml Sodium carbonate solution 30 g Glycine	
Mindesthaltbarkeit minimum shelf life		22.06.2016
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 08.10.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.