


**Certificate of Analysis**

Product Name: Quick Coomassie Stain  
Lot Number: 2025-004  
Appearance: Emerald Green Solution; May Be Slightly Turbid  
Absorbance @ 450nm (0.250-0.450): 0.350  
pH (1.00-1.20): 1.00

**Performance Analysis:**

1. 2µg Total Protein Load
2. Denaturing/Reducing Conditions on 4-12% Bis-Tris
3. BSA Stained With 60 Minutes Exposure with Quick Coomassie Stain.



Manufactured in UK under quality system management.

**George Storrs**  
Digitally signed  
by George Storrs  
Date: 2025.02.27  
15:15:25 Z

Signature

George Storrs  
Print name

Production Scientist  
Title

**Karl Culley**  
Digitally signed by Karl  
Culley  
Date: 2025.02.28  
10:52:55 Z

Signed by QA (or designee)

Calibre Scientific  
Unit 5a, R-Evolution, The Advanced Manufacturing Park, Rotherham, S60 5XA  
+44 (0) 114-224-2257  
info@proteinark.com | www.proteinark.com

Printed copies are "Uncontrolled Reference Only"  
It is the responsibility of the user to ensure that the published version is used.  
Confidential and Proprietary Information.