


Certificate of Analysis

Product Name: Quick Coomassie Stain
Lot Number: 2025-030
Appearance: Emerald Green Solution; May Be Slightly Turbid
Absorbance @ 450nm (0.250-0.450): 0.348
pH (1.00-1.20): 1.00

Performance Analysis:

- 2µg Total Protein Load
- Denaturing/Reducing Conditions on 4-12% Bis-Tris
- BSA Stained With 60 Minutes Exposure with Quick Coomassie Stain.



Manufactured in UK under quality system management.

George Storrs Digitally signed
by George Storrs
Date: 2025.09.12
08:35:33 +01'00'

Signature

George Storrs
Print name

Production Scientist
Title

Karl Culley Digitally signed by Karl Culley
Date: 2025.09.12 08:49:45
+01'00'

Signed by QA (or designee)

Calibre Scientific
Unit 5a, R-Evolution, The Advanced Manufacturing Park, Rotherham, S60 5XA
+44 (0) 114-224-2257
info@proteinark.com | www.proteinark.com

Printed copies are "Uncontrolled Reference Only"
It is the responsibility of the user to ensure that the published version is used.
Confidential and Proprietary Information.