


Certificate of Analysis

Product Name: Quick Coomassie Stain
Lot Number: 2025-035
Appearance: Emerald Green Solution; May Be Slightly Turbid
Absorbance @ 450nm (0.250-0.450): 0.355
pH (1.00-1.20): 1.00

Performance Analysis:

1. 2µg Total Protein Load
2. Denaturing/Reducing Conditions on 4-12% Bis-Tris
3. BSA Stained With 60 Minutes Exposure with Quick Coomassie Stain.



Manufactured in UK under quality system management.

Signature

George Storrs
Print name

Production Scientist
Title

Signed by QA (or designee)

Calibre Scientific
Unit 5a, R-Evolution, The Advanced Manufacturing Park, Rotherham, S60 5XA
+44 (0) 114-224-2257
info@proteinark.com | www.proteinark.com

Printed copies are "Uncontrolled Reference Only"
It is the responsibility of the user to ensure that the published version is used.
Confidential and Proprietary Information.

Distributed by SERVA Electrophoresis GmbH
SERVA Cat. no. 35081 Lot 251198