

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Siliconöl DC 200; 10 cst pract.	Kat.Nr./Cat.No. : 35132
Silicone DC 200 fluid; 10 cst, pract.	Lot/Contr.No. : 160582

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
Identität identity	IR	entspricht corresponds	entspricht corresponds
Dichte (g/ml) specific gravity	1013 hPa/ 25 °C	0.934 - 0.940	entspricht corresponds
Brechungsindex refractive index	25 °C	1.3980 - 1.4000	1.3992
Viskosität (cSt) viscosity		9.0 - 11.0	10.3
Flüchtige Bestandteile (%) volatile content		max. 1.0	0.13
Mindesthaltbarkeit minimum shelf life			17.11.2018
Lagerung (°C) storage			+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 09.06.2016

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.