

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

Siliconöl DC 200; 50 cst pract.	Kat.Nr./Cat.No. : 35134
Silicon DC 200 fluid; 50 cst pract.	Lot/Contr.No. : 251146

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, viskose Flüssigkeit clear, viscous liquid	entspricht corresponds
Identität identity	IR	entspricht corresponds	entspricht corresponds
Dichte (g/ml) specific gravity		0.957 – 0.963	entspricht corresponds
Brechungsindex refractive index	25 °C	1.4010 – 1.4030	1.4022
Viskosität (mm²/s) viscosity		47.5 – 52.5	51.4
Flüchtige Bestandteile (%) volatile content		max. 0.5	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			04.04.2028
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 08.01.2026

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.