

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Siliconöl DC 200; 350 cst pract.	Kat.Nr./Cat.No. : 35136
Silicon DC 200 fluid; 350 cst pract.	Lot/Contr.No. : 111277

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, viskose Flüssigkeit clear, viscous liquid	entspricht corresponds
Dichte (g/ml) specific gravity		0.966 – 0.972	entspricht corresponds
Brechungsindex refractive index	25 °C	1.4025 – 1.4045	1.4032
Viskosität (mm²/s) viscosity		332.5 – 367.5	347.9
Flüchtige Bestandteile (%) volatile content		max. 0.5	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			05/ 2015
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.08.2014

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.