

**A N A L Y S E N Z E R T I F I K A T**  
CERTIFICATE OF ANALYSIS

<b>Siliconöl DC 710; 500 cst pract.</b>	<b>Kat.Nr./Cat.No. : 35149</b>
Silicone DC 710 fluid; 500 cst pract.	<b>Lot/Contr.No. : 170312</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		klare, viskose Flüssigkeit clear, colorless liquid	entspricht corresponds
<b>Dichte (g/ml)</b> specific gravity	25 °C	1.10 – 1.11	entspricht corresponds
<b>Brechungsindex</b> refractive index	25 °C	1.5300 – 1.5380	1.5357
<b>Viskosität (cst)</b> viscosity	25 °C	450.0 – 550.0	536.9
<b>Mindesthaltbarkeit</b> minimum shelf life			20.10.2021
<b>Lagerung (°C)</b> storage			+15 bis / to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 30.03.2017**

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.