

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Siliconöl DC 710; 500 cst pract.	Kat.Nr./Cat.No. : 35149
Silicone DC 710 fluid; 500 cst pract.	Lot/Contr.No. : 180388

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, viskose Flüssigkeit clear, colorless liquid	entspricht corresponds
Dichte (g/ml) specific gravity	25 °C	1.10 – 1.11	entspricht corresponds
Brechungsindex refractive index	25 °C	1.5300 – 1.5380	1.5342
Viskosität (cst) viscosity	25 °C	450.0 – 550.0	452
Mindesthaltbarkeit minimum shelf life			24.01.2022
Lagerung (°C) storage			+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 12.04.2018

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.