

CERTIFICATE OF ANALYSIS

Thrombin from bovine plasma min. 1000 units/mg protein lyophil. E.C. 3.4.21.5	Cat. No. : 36402 Lot. No. : 150112
---	---

Parameter	Method	Specification	Result
Molecular weight		ca. 37 000	
Appearance		white lyophilisate	corresponds
Specific Activity (NIH U/mg protein)	Thrombin	min. 1000	2402
Activity (U/mg solid)	Thrombin		249.8
Minimum shelf life			01/ 2018
Storage (°C)			+2 to +8

Unit definition:

One NIH unit clots a standard fibrinogen solution in 15 sec. at 37°C.

This material is derived from animals born and sourced from the United States. The animals are from a natural beef program in accordance to FDA Regulation 21 CFR 589.2000, which prohibits the use of animal proteins in ruminant feed.

The animals, both pre- and post-mortem, were apparently free from infectious and contagious diseases. Therefore to the best of our knowledge and belief, we can state that BSE/TSE should not be a concern.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Quality Control

Printing date: 28.01.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.