

CERTIFICATE OF ANALYSIS

Thrombin from bovine plasma	Cat. No. : 36402
min. 1000 units/mg protein lyophil. E.C. 3.4.21.5	Lot. No. : 250128

Parameter	Method	Specification	Result
Molecular weight		ca. 37 000	
Appearance		white lyophilisate	corresponds
Specific Activity (NIH U/mg protein)	Thrombin	min. 1000	2750
Activity (U/mg solid)	Thrombin		895
Origin			New Zealand
Minimum shelf life			11.02.2028
Storage (°C)			+2 to +8

Unit definition:

One NIH unit clots a standard fibrinogen solution in 15 sec. at 37 °C.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Quality Control

Printing date: 13.02.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.