

CERTIFICATE OF ANALYSIS



PNGase F, recombinant	Cat.No.: : 36404
	Contr.No.: 250192

Parameter	Method	Specification	Result
Molecular weight	PNGase F	ca. 36 000	
Purity (%)	SDS-PAGE	≥ 95	≥ 95
Activity (U/ml)		1 000 000	1 000 000
Form	The enzyme is supplied as a liquid in 1 x PBS (Phosphate Buffered Saline; 137 mM NaCl, 10 mM Phosphate, 2.7 mM KCl, pH 7.4)		
Origin		E. coli	corresponds
Minimum shelf life			27.02.2027
Storage (°C)			-15 bis/ to -25

Unit definition:

Achieves complete deglycosylation of 10 µg of RNase B incubated in 1 x PBS with 1 µl of PNGase F for 5 – 10 minutes at 37 °C or room temperature. Separation of reaction products are visualized by SDS-PAGE.

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 28.02.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.