

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>Tris(hydroxymethyl)aminomethan</b>	<b>Kat.Nr./Cat.No. :</b>	<b>37181</b>
TRIS electrophoresis grade	<b>Lot/Contr.No. :</b>	<b>141187</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Summenformel / formula</b>		C <sub>4</sub> H <sub>11</sub> NO <sub>3</sub>	
<b>Molekulargewicht</b> molecular weight		121.1	
<b>Aussehen</b> appearance		weiße Kristalle white crystals	entspricht corresponds
<b>Gehalt (%)</b> assay	Titration	min. 99.9	100.1
<b>Identität</b> identity	IR	entspricht corresponds	entspricht corresponds
<b>Trocknungsverlust (%)</b> loss on drying	3 h, 105 °C	max. 1.0	0.03
<b>pH</b>	5 % in H <sub>2</sub> O	10.0 – 11.5	10.7
<b>Löslichkeit</b> solubility	10 % in H <sub>2</sub> O	klar und farblos clear and colorless	entspricht corresponds
<b>A 1 cm/ 10 % in H<sub>2</sub>O</b>	235 nm	max. 0.60	0.104
	260 nm	max. 0.15	0.020
	280 nm	max. 0.15	0.022
	430 nm	max. 0.07	0.002
<b>Schmelzpunkt (°C)</b> melting point		168 – 172	171
<b>Glührückstand (%)</b> residue on ignition		max. 0.1	< 0.1
<b>Schwermetalle (ppm)</b> heavy metals	als/ as Pb	max. 1	< 1
<b>Eisen (ppm)/ iron</b>	Fe	max. 1	< 1
<b>Anwendung</b> application	Gelelektrophorese gel electrophoresis	entspricht corresponds	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			01/ 2019
<b>Lagerung (°C)/ storage</b>			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 19.11.2014**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.