

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

N-Tris(hydroxymethyl)methyl-3-aminopropan-sulfonsäure p. A.	Kat.Nr./Cat.No. :	37194
TAPS analytical grade	Lot/Contr.No. :	060853

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₇ H ₁₇ NO ₆ S	
Molekulargewicht molecular weight		243.28	
Aussehen appearance		weißes Pulver white powder	entspricht corresponds
Gehalt (%) assay	Titration	min. 99.0	99.9
Schmelzpunkt (°C) melting point		238 – 245 (Zersetzung/ decomp.)	238 – 239
A 1 cm/ 230 nm	10 % in H ₂ O	max. 0.2	0.035
A 1 cm/ 260 nm	10 % in H ₂ O	max. 0.07	0.022
A 1 cm/ 280 nm	10 % in H ₂ O	max. 0.07	0.019
pH	10 % in H ₂ O	4.5 – 6.5	5.2
Identität identity	IR	entspricht corresponds	entspricht corresponds
Schwermetalle (ppm) heavy metals	als/ as Pb	max. 10	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			12/ 2012
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 15.12.2011

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.