

CERTIFICATE OF ANALYSIS

<b>Trypsin 1:250 from bovine pancreas</b> lyophil.	<b>Cat.No. :</b> 37289
	<b>Contr.No. :</b> 190039

Parameter	Method	Specification	Result
<b>Appearance</b>		white lyophilisate	corresponds
<b>Activities:</b>			
<b>Trypsin (U/mg)</b>	Bergmeyer		4.2
<b>Trypsin (NF/USP U/mg)</b>	USP	≥ 250	349
<b>Chymotrypsin (NF/USP U/mg)</b>	USP	≤ 75	59
<b>Microbiological analysis</b>	Total viable count	≤ 10 000 cfu/g	corresponds
	Pseudomonas aeruginosa	negative	corresponds
	Salmonella species	negative	corresponds
	Staphylococcus aureus	negative	corresponds
	E.coli	negative	corresponds
<b>Application</b>		suitable for cell culture	
<b>Minimum shelf life</b>			29.11.2020
<b>Storage (°C)</b>			-15 bis/ to -25

**Unit definition**

**Trypsin (Bergmeyer)**

1 U catalyzes the hydrolysis of 1 µmol N-α-benzoyl-L-arginine ethyl ester (BAEE) per minute at 25°C, pH 8.0.

**We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.**

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 17.01.2019

Christian Monsler

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.