

CERTIFICATE OF ANALYSIS

Trypsin inhibitor from soybean lyophil.	Cat.No. : 37329 Contr.No.: 160375
---	--

Parameter	Method	Specification	Result
Molecular weight		ca. 22 000	
Appearance		pale yellow lyophilisate	corresponds
Solubility	1 mg/ml in 0.067 M potassium-phosphate buffer, pH 7.6	clear and colorless	corresponds
Activity (U/mg)		min. 13 000	15 000
Minimum shelf life			08.08.2019
Storage (°C)			+2 to +8

Unit definition

1 IU (inhibitor units) inhibits 1 U trypsin as defined by cleavage of 1 mmole BAEE (N-benzoyl-L-arginine ethyl ester) per minute.

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 23.07.2018

Christian Monsler

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.