

CERTIFICATE OF ANALYSIS

Urease from Jack bean lyophil.	Cat.No. : 37799
	Contr.No.: 140028

Parameter	Method	Specification	Result
Molecular weight		ca. 545 000	
Appearance		white lyophilisate	corresponds
Activity (U/mg)	Urease by Nessler's	min. 220	312
Activity in other units (U/mg)	Urease by Bergmeyer	min. 90	130
Minimum shelf life			01/ 2016
Storage (°C)			+2 to +8

Unit definitions

Urease (by Nessler's)

1 Unit is that amount of enzyme causing the liberation of 1 µmol ammonia per minute at 25°C, pH 7.0.

Urease (by Bergmeyer)

1 Unit is that amount of enzyme causing the hydrolysis of 1 µmol of urea per minute at 25°C, pH 8.0.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Quality Control

Printing date: 13.01.2014

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.