

CERTIFICATE OF ANALYSIS

Protein Molecular Weight Standards	Cat. No. :	39064
range 12.4 – 450 kDa for native PAGE	Lot. No. :	P160064

Protein Components	Molecular weight
Cytochrome C horse heart	12 400
Myoglobin equine	17 800
Chymotrypsinogen A	25 000
Albumin egg	45 000
Albumin bovine (BSA)	67 000
Aldolase rabbit	160 000
Catalase bovine	240 000
Ferritin horse	450 000
Storage (°C)	-15 to -25
Minimum shelf life	13.12.2017

The kit MS II contains 25 mg of each of the above mentioned, water soluble substances with well defined molecular weights in native electrophoresis.

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 28.11.2016

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

This report has been computer generated and does not contain a signature.