

# CERTIFICATE OF ANALYSIS

<b>Protein Molecular Weight Standards</b>	<b>Cat. No. :</b>	<b>39064</b>
range 12.4 – 450 kDa for native PAGE	<b>Lot. No. :</b>	<b>P170317</b>

<b>Protein Components</b>	<b>Molecular weight</b>
Cytochrome C horse heart	12 400
Myoglobin equine	17 800
Chymotrypsinogen A	25 000
Albumin egg	45 000
Albumin bovine (BSA)	67 000
Aldolase rabbit	160 000
Catalase bovine	240 000
Ferritin horse	450 000
Storage (°C)	-15 to -25
Minimum shelf life	30.07.2019

The kit MS II contains 25 mg of each of the above mentioned, water soluble substances with well defined molecular weights in native electrophoresis.

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date: 25.05.2018**

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

This report has been computer generated and does not contain a signature.