

CERTIFICATE OF ANALYSIS



Protein Molecular Weight Standards range 12.3 – 450 kDa for native PAGE	Cat. No. : 39064 Lot. No. : P240234
---	--

Protein Components	Molecular weight
Cytochrome C (porcine, heart)	12 300
Myoglobin (equine)	17 800
Chymotrypsinogen A	25 000
Albumin (egg)	45 000
Albumin (bovine, BSA)	67 000
Aldolase (rabbit)	160 000
Catalase (bovine)	240 000
Ferritin (horse)	450 000
Storage (°C)	-15 to -25
Minimum shelf life	02.10.2025

The kit contains 25 mg of each of the above mentioned, water soluble substances with well defined molecular weights in native electrophoresis.

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 21.10.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.