

CERTIFICATE OF ANALYSIS**Protease-Inhibitor-Mix P****Cat.No. : 39103****Lot No.: P120027**

Cat. number	Components	Target Protease
12745	AEBSF-HCl	Serine Proteases
21100	E-64	Cysteine Proteases
52682	Pepstatin A	Aspartic Proteases
32010	1,10-Phenanthroline	Metalloproteases
14980	Bestatin	Aminopeptidase B and Leucine Aminopeptidase
51867	Leupeptin	Cysteine Proteases and Trypsin-like Serine Proteases
Storage (°C)		-15 to -25
Minimum shelf life		01/ 2015

Reconstitution

Dissolve the content of the vial (64.2 mg) in 1ml DMSO to obtain a 100x stock solution (suitable for the treatment of 100 ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for at least 4 weeks).

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control**Printing date:** 19.01.2012

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

This report has been computer generated and does not contain a signature.