

CERTIFICATE OF ANALYSIS

Protease-Inhibitor-Mix P

Cat.No. : 39103

Lot No.: P240054

Cat. number	Components	Target Protease
12745	AEBSF-HCl	Serine Proteases
21100	E-64	Cysteine Proteases
52682	Pepstatin A	Aspartic Proteases
32010	1,10-Phenanthroline	Metalloproteases
14980	Bestatin	Aminopeptidase B and Leucine Aminopeptidase
51867	Leupeptin	Cysteine Proteases and Trypsin-like Serine Proteases
Storage (°C)		-15 to -25
Minimum shelf life		19.02.2027

Reconstitution

Dissolve the content of the vial (64.2 mg) in 1ml DMSO to obtain a 100x stock solution (suitable for the treatment of 100 ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for at least 4 weeks).

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 23.02.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.