

CERTIFICATE OF ANALYSIS

Protease-Inhibitor-Mix FY

Cat.No. : 39104

Lot No. : P160138

| Cat. number | Components | Target Protease |
|--------------------|---------------------|------------------------|
| 12745 | AEBSF-HCl | Serine Proteases |
| 21100 | E-64 | Cysteine Proteases |
| 52682 | Pepstatin A | Aspartic Proteases |
| 32010 | 1,10-Phenanthroline | Metalloproteases |

Storage (°C) -15 to -25

Minimum shelf life 23.03.2019

Reconstitution Dissolve the content of the vial (62.6 mg) in 1 ml DMSO to obtain a 100x stock solution (suitable for the treatment of 100 ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for at least 4 weeks).

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 24.03.2016

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.