

CERTIFICATE OF ANALYSIS

Protease-Inhibitor-Mix FY

Cat.No. : 39104

Lot No. : P170212

Cat. number	Components	Target Protease
12745	AEBSF-HCl	Serine Proteases
21100	E-64	Cysteine Proteases
52682	Pepstatin A	Aspartic Proteases
32010	1,10-Phenanthroline	Metalloproteases

Storage (°C) -15 to -25

Minimum shelf life 14.06.2020

Reconstitution Dissolve the content of the vial (62.6 mg) in 1 ml DMSO to obtain a 100x stock solution (suitable for the treatment of 100 ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for at least 4 weeks).

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 22.06.2017

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

This report has been computer generated and does not contain a signature.