

**CERTIFICATE OF ANALYSIS****Protease-Inhibitor-Mix B****Cat.No. : 39105****Lot No.: P100384**

<b>Cat. number</b>	<b>Components</b>	<b>Target Protease</b>
12745	AEBSF-HCl	Serine Proteases
21100	E-64	Cysteine Proteases
52682	Pepstatin A	Aspartic Proteases
14980	Bestatin	Aminopeptidase B and Leucine Aminopeptidase
11280	EDTA-Disodium	Metalloproteases
<b>Minimum shelf life</b>		11/ 2014
<b>Storage (°C)</b>		-15 to -25
<b>Reconstitution</b>	Dissolve the content of the vial (26.4mg) in 1ml 90% DMSO to obtain a 100x stock solution (suitable for the treatment of 100ml tissue extract). Slight turbidity in the reconstituted solution is normal. After reconstitution, aliquot and freeze at -20°C (stable for at least 4 weeks).	

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control****Printing date:** 16.09.2013

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.