

CERTIFICATE OF ANALYSIS

Protease-Inhibitor-Mix B

Cat.No. : 39105

Lot No.: P110162

Cat. number	Components	Target Protease
12745	AEBSF-HCl	Serine Proteases
21100	E-64	Cysteine Proteases
52682	Pepstatin A	Aspartic Proteases
14980	Bestatin	Aminopeptidase B and Leucine Aminopeptidase
11280	EDTA-Disodium	Metalloproteases
Minimum shelf life		05/ 2015
Storage (°C)		-15 to -25
Reconstitution	Dissolve the content of the vial (26.4 mg) in 1 ml 90 % DMSO to obtain a 100x stock solution (suitable for the treatment of 100 ml tissue extract). Slight turbidity in the reconstituted solution is normal. After reconstitution, aliquot and freeze at -20 °C (stable for at least 4 weeks).	

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 25.03.2014

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.