

CERTIFICATE OF ANALYSIS

Protease-Inhibitor-Mix B

Cat.No. : 39105

Lot No.: P210409

| Cat. number | Components | Target Protease |
|---------------------------|--|--|
| 12745 | AEBSF-HCl | Serine Proteases |
| 21100 | E-64 | Cysteine Proteases |
| 52682 | Pepstatin A | Aspartic Proteases |
| 14980 | Bestatin | Aminopeptidase B and Leucine Aminopeptidase |
| 11280 | EDTA-Disodium | Metalloproteases |
| Minimum shelf life | | 18.11.2025 |
| Storage (°C) | | -15 to -25 |
| Reconstitution | Dissolve the content of the vial (26.4 mg) in 1 ml 90 % DMSO to obtain a 100x stock solution (suitable for the treatment of 100 ml tissue extract). Slight turbidity in the reconstituted solution is normal. After reconstitution, aliquot and freeze at -20 °C (stable for at least 4 weeks). | |

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 10.10.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.