

# CERTIFICATE OF ANALYSIS

**Protease-Inhibitor-Mix B**

**Cat.No. : 39105**

**Lot No.: P240214**

<b>Cat. number</b>	<b>Components</b>	<b>Target Protease</b>
12745	AEBSE-HCl	Serine Proteases
21100	E-64	Cysteine Proteases
52682	Pepstatin A	Aspartic Proteases
14980	Bestatin	Aminopeptidase B and Leucine Aminopeptidase
11280	EDTA-Disodium	Metalloproteases
<b>Minimum shelf life</b>		02.07.2027
<b>Storage (°C)</b>		-15 to -25
<b>Reconstitution</b>	Dissolve the content of the vial (26.4 mg) in 1 ml 90 % DMSO to obtain a 100x stock solution (suitable for the treatment of 100 ml tissue extract). Slight turbidity in the reconstituted solution is normal. After reconstitution, aliquot and freeze at -20 °C (stable for at least 4 weeks).	

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**

**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 18.07.2024**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**

This report has been computer generated and does not contain a signature.