

# Certificate of Analysis

**SERVA Native Marker**  
Liquid Mix

**Cat.No. : 39213**  
**Lot No.: P250211**

<b>Protein Components</b>	<b>Molecular weight (kDa)</b>
<b>Ferritin</b> (horse)	720 / 450
<b>Phycoerythrin B</b>	240
<b>Lactate dehydrogenase</b> (porcine heart)	146
<b>Albumin</b> (bovine)	67
<b>Albumin</b> (egg)	45
<b>Trypsin inhibitor</b> (soybean)	21
<b>Storage</b> (°C)	-15 to -25
<b>Minimum shelf life</b>	26.06.2026

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 03.07.2025

Dr. Judith Koch

Jan Dennis Lorenzen

This report has been computer generated and does not contain a signature.