

CERTIFICATE OF ANALYSIS



SERVA Unstained Protein Standard 6.5 - 97 kDa

Cat.No. : 39214

Lot No. : P240317

Protein Components	Molecular weight
Phosphorylase B	97 400
Albumin bovine (BSA)	67 000
Albumin egg	45 000
Carbonic anhydrase	29 000
Trypsin inhibitor (soybean)	21 000
Cytochrome C	12 500
Trypsin inhibitor (bovine lung)	6 500
Minimum shelf life	15.10.2025
Storage (°C)	-15 to -25

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

**SERVA Electrophoresis GmbH
Quality Control**

Printing date: 17.10.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.