

Certificate of Analysis

SERVA Native Marker	Cat.No. :	39219
Liquid Mix	Lot No.:	P150102

Protein Components	Molecular weight (kDa)
Ferritin (horse)	720 / 450
Urease (jack bean)	545 / 272
Lactate dehydrogenase (porcine heart)	146
Albumin (bovine)	67
Albumin (egg)	45
Trypsin inhibitor (soybean)	21
Storage (°C)	-15 to -25
Minimum shelf life	12.06.2016

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 17.11.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.