

Certificate of Analysis

SERVA Native Marker	Cat.No. :	39219
Liquid Mix	Lot No.:	P180121

Protein Components	Molecular weight (kDa)
Ferritin (horse)	720 / 450
Urease (jack bean)	545 / 272
Lactate dehydrogenase (porcine heart)	146
Albumin (bovine)	67
Albumin (egg)	45
Trypsin inhibitor (soybean)	21
Storage (°C)	-15 to -25
Minimum shelf life	03.07.2019

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 06.04.2018

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

This report has been computer generated and does not contain a signature.