

CERTIFICATE OF ANALYSIS

SingleQuant Assay Kit	Cat.No. : 39226
	Contr.No. : P110499

Kit Components	Storage
Reference standard (BSA)	aliquots at -20°C
Wash solution	room temperature
Elution solution	room temperature
SingleQuant stock solution	+2 to +8°C
Application	protein determination
Minimum shelf life	11/ 2013

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 29.11.2011

Dipl.-Ing. (FH) Bernhard Göckel

Dipl.-Ing. (FH) Christian Nowak

This report has been computer generated and does not contain a signature.