

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

BCA Protein Assay Micro Kit	Kat.Nr./Cat.No. : 39229
BCA Protein Assay Micro Kit	Lot/Contr.No. : P250134

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Anwendung application	BCA Protein Assay	entspricht corresponds	entspricht corresponds
Komponenten components		Reagent A: Puffer-Lösung Reagent B: BCA-Lösung Reagent C: Kupfersulfat-Lösung Proteinstandard (Rinderserumalbumin)	entspricht
		Reagent A: Buffer-Lösung Reagent B: BCA solution Reagent C: copper sulphate solution Protein standard (bovine serum albumin)	corresponds
Mindesthaltbarkeit minimum shelf life			04.08.2026
Lagerung (°C) storage			+2 bis/ to +8

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 10.04.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.