

**A N A L Y S E N Z E R T I F I K A T**  
CERTIFICATE OF ANALYSIS

<b>SERVA Purple Protein Quantifizierungs-Assay</b>	<b>Kat.Nr./Cat.No. :</b>	<b>39235</b>
SERVA Purple Protein Quantification Assay	<b>Lot/Contr.No. :</b>	<b>P190157</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		gelbe Flüssigkeit yellow liquid	entspricht corresponds
<b>Reinheit (%)</b> purity	HPLC	Status	> 95
<b>Mindesthaltbarkeit</b> minimum shelf life			16.04.2021
<b>Lagerung (°C)</b> storage			-15 bis/ to -25

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 18.04.2019

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.