

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

SERVA Purple Protein Quantifizierungs-Assay	Kat.Nr./Cat.No. :	39235
SERVA Purple Protein Quantification Assay	Lot/Contr.No. :	P190157

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		gelbe Flüssigkeit yellow liquid	entspricht corresponds
Reinheit (%) purity	HPLC	Status	> 95
Mindesthaltbarkeit minimum shelf life			16.04.2022
Lagerung (°C) storage			-15 bis/ to -25

SERVA Electrophoresis GmbH

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.05.2021

Christian Monsler

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.