

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

RIPA-Puffer	Kat.Nr./Cat.No. : 39244
RIPA Buffer	Lot/Contr.No. : P180391

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	klare, farblose Lösung clear, colorless solution	entspricht corresponds
Zusammensetzung composition		
Tris buffer pH 7.5, 1 M solution	25 mM	entspricht/ corresponds
NaCl	150 mM	entspricht/ corresponds
Triton® X-100	1 %	entspricht/ corresponds
Deoxycholic acid-Na-salt	1 %	entspricht/ corresponds
SDS solution, 20 %	0.1 %	entspricht/ corresponds
Mindesthaltbarkeit minimum shelf life		03.12.2020
Lagerung (°C) storage		+2 bis/ to +8

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 10.12.2018

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.