

# CERTIFICATE OF ANALYSIS

---

**SERVA Unstained Protein Standard IV**

**Cat.No. : 39250**

**Lot No. : 200730**

---

<b>Protein Recombinant</b>	<b>Molecular Weight (kDa)</b>
<b>Protein 1</b>	200
<b>Protein 2</b>	150
<b>Protein 3</b>	100
<b>Protein 4</b>	85
<b>Protein 5</b>	60
<b>Protein 6</b>	50
<b>Protein 7</b>	40
<b>Protein 8</b>	30
<b>Protein 9</b>	25
<b>Protein 10</b>	20
<b>Protein 11</b>	15
<b>Protein 12</b>	10
<b>Storage (°C)</b>	-15 to -25
<b>Minimum Shelf Life</b>	20.05.2022

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 22.07.2020**

Christian Monsler

Marina Konter

This report has been computer generated and does not contain a signature.