

CERTIFICATE OF ANALYSIS

SERVA Unstained Protein Standard IV

Cat.No. : 39250

Lot No. : 210165

Protein Recombinant	Molecular Weight (kDa)
Protein 1	200
Protein 2	150
Protein 3	100
Protein 4	85
Protein 5	60
Protein 6	50
Protein 7	40
Protein 8	30
Protein 9	25
Protein 10	20
Protein 11	15
Protein 12	10
Storage (°C)	-15 to -25
Minimum Shelf Life	15.02.2023

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 18.02.2021

Christian Monsler

Patricia Milford

This report has been computer generated and does not contain a signature.