

CERTIFICATE OF ANALYSIS

SERVA Triple Color Protein Standard I
15 – 180 kDa, Liquid Mix

Cat.No. : 39251
Lot No. : 211183

Protein Components	Molecular Weight (kDa)
Protein 1 (blue)	180
Protein 2 (blue)	130
Protein 3 (blue)	100
Protein 4 (red)	70
Protein 5 (green)	60
Protein 6 (blue)	40
Protein 7 (red)	30
Protein 8 (blue)	20
Protein 9 (green)	15
Storage (°C)	-15 to -25
Minimum Shelf Life	06.01.2024

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 13.10.2021

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.