

CERTIFICATE OF ANALYSIS



SERVA Triple Color Protein Standard I
15 – 170 kDa, Liquid Mix

Cat.No. : 39251
Lot No. : 240707

Protein Components	Molecular Weight (kDa)
Protein 1 (blue)	170
Protein 2 (blue)	130
Protein 3 (blue)	100
Protein 4 (red)	70
Protein 5 (green)	60
Protein 6 (blue)	45
Protein 7 (blue)	35
Protein 8 (blue)	25
Protein 9 (green)	15
Storage (°C)	-15 to -25
Minimum Shelf Life	24.01.2027

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 07.08.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.