

**CERTIFICATE OF ANALYSIS****SERVA Triple Color Protein Standard II**  
11 – 180 kDa, Liquid Mix**Cat.No. : 39257**  
**Lot No. : 130928**

<b>Protein Components</b>	<b>Molecular Weight (kDa)</b>
<b>Protein 1 (blue)</b>	180
<b>Protein 2 (blue)</b>	135
<b>Protein 3 (blue)</b>	100
<b>Protein 4 (orange)</b>	75
<b>Protein 5 (blue)</b>	63
<b>Protein 6 (blue)</b>	48
<b>Protein 7 (blue)</b>	35
<b>Protein 8 (green)</b>	25
<b>Protein 9 (blue)</b>	17
<b>Protein 10 (blue)</b>	11
<b>Storage (°C)</b>	-15 to -25
<b>Minimum Shelf Life</b>	08/ 2015

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control****Druckdatum/ printing date: 16.08.2013**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Steinmann

This report has been computer generated and does not contain a signature.